
BOSTON HEART DIAGNOSTICS Beta Cell Function and Risk Index


The **Beta Cell Function and Risk Index** is a calculation designed to give insights into beta cell function. It allows clinicians to assess baseline beta cell function, consider treatment options, and follow changes over time.

IMPROVING BETA CELL FUNCTION CAN PREVENT OR DELAY TYPE 2 DIABETES

PROGRESSIVE BETA CELL FAILURE

- In patients with chronic insulin resistance (low HOMA-S), beta cell activity can be very high (increased HOMA-B) in order to stabilize glucose levels.
- Over time, beta cells can become progressively impaired, potentially resulting in prediabetes followed by type 2 diabetes. While HOMA-S levels remain low, HOMA-B levels gradually go from high to low as beta cell activity deteriorates.
- If HOMA-B levels are very low, then supplemental insulin may be needed due to loss of beta cell function.
- There are multiple treatment strategies for slowing or reversing progressive beta cell dysfunction, prediabetes, and type 2 diabetes.⁴⁻⁸






ABOUT BETA CELL FUNCTION AND RISK INDEX

- Beta cells in the pancreas secrete insulin to regulate blood glucose levels.
- The Homeostasis Model Assessment can be used to estimate Beta Cell Function (HOMA-B), Insulin Sensitivity (HOMA-S) and Insulin Resistance (HOMA-IR). These are calculations based on the relationship between fasting insulin and glucose.¹⁻³
- Clinical decision-making can be facilitated by characterizing a patient's beta cell activity after accounting for the degree of insulin sensitivity.
- The Boston Heart Beta Cell Function & Risk Index is a tool for patient engagement and motivation



QUANTIFYING BETA CELL FUNCTION AND RISK



For any specific level of insulin sensitivity, patients with known diabetes tend to have relatively weak beta cell response compared to those with normal glucose metabolism. Patients with prediabetes fall in between.

For example, a HOMA-B of 100% would be normal for an individual with normal insulin sensitivity (HOMA-S = 100%), but would be abnormally low for an individual with low insulin sensitivity (HOMA-S = 25%), consistent with the point where prediabetes progresses to type 2 diabetes.

ORDERING AND REPORTING


Ordering Information

• The order code is 408C (includes Insulin, Glucose, HOMA-IR, HOMA-S, HOMA-B, Beta Cell Function and Risk Score).

Reporting

- The Beta Cell Function and Risk Index will be reported to the right of the Metabolic section of the lab report.
- To facilitate clinical interpretation, the Boston Heart Laboratory report provides the HOMA-IR, HOMA-S, HOMA-B values, and a sliding scale with a percentile risk index value between 0 and 100.
- Interpretations provide commentary regarding the degree of metabolic risk, insulin sensitivity, and beta cell activity. Treatment considerations include lifestyle modification (weight loss, exercise, limitation of refined carbohydrates), and evidence-based medication options when indicated.
- Follow-up test results indicate the degree and direction of change since the last assessment.

bostonheartdiagnostics.com



Boston Heart Beta Cell Function and Risk Index
Optimal Borderline Increased Risk
Current
Previous

Interpretation: Consider prediabetes and metabolic risk with low insulin sensitivity and very elevated beta cell activity. Beta cell risk is in the 50th percentile. The risk has moderately decreased from the previous result (60th Percentile).

Consideration: Consider weight loss, exercise, and/or metformin to increase insulin sensitivity. For elevated beta cell activity (HOMA-B), consider limiting refined carbohydrates such as sugars and starches to reduce insulin demand.

References

877.425.1252

- 1. Wallace TM, et al. Use and abuse of HOMA modeling. Diabetes Care 2004 Jun;27(6):1487-95.
- Hill NR, et al. Expansion of the homeostasis model assessment of -cell function and insulin resistance to enable clinical trial outcome modeling through the interactive adjustment of physiology and treatment effects: iHOMA2.Diabetes Care. 2013 Aug;36(8):2324-30.

increased risk of developing diabetes, hyperlipidemia, hypertension and CVD. Based on the HbA1c value, the estimated Average Glucose (eAG) is 120 mg/dL which includes the non-fasting state.

Consideration: Consider encouraging dietary modification supported by education and consider glucose lowering and/or insulin sensitizing medications. If indicated encourage weight reduction, smoking cessation, increased activity and control blood pressure.

- 3. Song Y, et al. Insulin sensitivity and insulin secretion determined by homeostasis model assessment and risk of diabetes in a multiethnic cohort of women: the Women's Health Initiative Observational Study. Diabetes Care. 2007 Jul;30(7):1747-52
- 4. Sheng Z, et al. Effects of lifestyle modification and anti-diabetic medicine on prediabetes progress: A systematic review and meta-analysis. Front Endocrinol (Lausanne). 2019;10:455.
- 5. Hu Y, et al. Short-term intensive therapy in newly diagnosed type 2 diabetes partially restores both insulin sensitivity and -cell function in subjects with long-term remission Diabetes Care. 2011
- Aug;34(8):1848-53
- 6. Wang H, et al. Predictors of long-term glycemic remission after 2-week intensive insulin treatment in newly diagnosed type 2 diabetes. J Clin Endocrinol Metab. 2019 Jun 1;104(6):2153-2162
- Lean MEJ, et al. Durability of a primary care-led weight-management intervention for remission of type 2 diabetes: 2-year results of the DiRECT open-label, cluster-randomised trial. Lancet Diabetes Endocrinol. 2019 May;7(5):344-355
- 8. Al-Mrabeh et al. Hepatic lipoprotein export and remission of human type 2 diabetes after weight loss. Cell Metab. 2020 Feb 4;31(2):233-249.e4
- 9. National Diabetes Statistics Report, 2020. Atlanta, GA: Centers for Disease Control and Prevention, US Department of Health and Human Services; 2020.



© 2020 Boston Heart Diagnostics Corporation. All rights reserved. The Boston Heart logo is a registered trademark of Boston Heart Diagnostics Corporation in the U.S. and in other countries. Boston Heart Diagnostics Corporation reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors. 03310320

