BOSTON HEART HDL MAP® ASSESSING AND TREATING SUBPOPULATION ABNORMALITIES

LARGE **a-1** HDL

Large α -mobility HDL particles are associated with healthy HDL maturation and decreased CVD risk. High levels of α -1 HDL are a marker of protection from heart disease.

Model of Spherical Lipoprotein developed by Boston Heart Diagnostics, composed of protein (yellow), phospholipids (blue), and free cholesterol (green) on the surface, and cholesteryl ester (dark green) and triglyceride (purple) in the core.



BOSTON HEART HDL MAP®

FORMATION OF HDL PARTICLES

HDL particles begin as free circulating apoA-I, which is made in the liver and intestine. Lipid free apoA-I gains phospholipids and becomes a very small HDL particle, called preB-1. These particles begin to collect free cholesterol and additional phospholipids from cells and progressively grow into larger HDL particles (small α -4, medium α -3, large α -2 and very large α -1). The conversion of free cholesterol on the HDL particle surface into cholesteryl ester results in the growth of HDL particles as the cholesteryl ester moves into the HDL particle core. When HDL particles become large α -2 or very large α -1 they dump cholesteryl ester to the liver or to triglyceride rich lipoproteins, the apoA-I then recycles back as a small HDL particle again.



BOSTON HEART HDL MAP®

OVERVIEW OF THE HDL MAP

Standard HDL-C tests only measure the total amount of HDL-C contained by all HDL particles. Boston Heart's exclusive HDL Map test measures the amount of apoA-I in the five most significant subpopulations, resulting in a deeper understanding of a patient's CVD risk.¹⁻⁴ This exclusive method provides an accurate and consistent indication of reverse cholesterol transport by separating the larger, cardioprotective particles (α -1) most associated with decreased CVD risk from the smaller HDL particles which have been associated with increased risk.¹⁻⁴

SCIENCE OF THE HDL MAP

The Boston Heart HDL Map test analyzes the distribution of HDL subpopulations in plasma using a proprietary gel electrophoresis technique which enables precise differentiation of HDL subparticles based on size.

Subsequent immunoblotting quantifies the amount of apoA-I, the main protein of HDL, in each of the five most important HDL subpopulations (very large α -1, large α -2, medium α -3, small α -4 and very small preB-1), providing more accurate disease characterization than HDL-C alone.²



BOSTON HEART HDL MAP[®]

TREATMENT ALGORITHM FOR HDL SUBPOPULATION ABNORMALITIES



ASSESS FOR ABNORMAL TEST RESULTS

MALE FEMALE	Low apoA-I in $\alpha\text{-1}$ HDL level <25 mg/dL Low apoA-I in $\alpha\text{-1}$ HDL level <35 mg/dL
MALE FEMALE	Low high-density lipoprotein cholesterol (HDL-C) level <40 mg/dL Low high-density lipoprotein cholesterol (HDL-C) level <50 mg/dL
ALL	High triglyceride (TG) level >150 mg/dL

ADDITIONAL CAUSES OF ELEVATED TRIGLYCERIDES⁶



ASSESS FOR ADDITIONAL CAUSES OF LIPID ABNORMALITIES

ADDITIONAL CAUSES OF LOW HDL CHOLESTEROL⁵

LIFESTYLE	CONDITIONS	MEDICATIONS	LIFESTYLE	CONDITIONS	MEDICATIONS
Cigarette smoking High sugar intake Excess caloric intake Excess alcohol intake (>2 drinks/day) High trans fat intake Physical inactivity	Common Causes - Elevated triglycerides [®] * - Diabetes mellitus* - Insulin resistance - Overweight/obesity - Kidney & liver dysfunction Other Causes - HIV - Polycystic ovarian syndrome - Acute or chronic inflammation - Hypothyroidism	 Non-selective beta blockers Androgenic steroids Progestins Isotretinoin Paradoxical response to fenofibrate and thiazolidinediones 	 Cigarette smoking High sugar intake Excess caloric intake Excess alcohol intake (>2 drinks/day) Physical inactivity High saturated fat intake, fried foods 	Common Causes Overweight/obesity Insulin resistance Diabetes mellitus Alcoholism Kidney dysfunction Other Causes HIV Cushing's disease Pregnancy Hypothyroidism 	 Non-selective beta blockers Thiazide diuretics Androgenic steroids Oral estrogens Oral contraceptives Isotretinoin Protease inhibitors Cyclosporin Glucocorticosteroid

*Condition must be optimally controlled to successfully increase HDL.



Assess for family history of an HDL disorder and premature heart disease

- Collect family history about premature cardiovascular disease (CVD) identified as male first degree relative < age 55 and female first degree relative < age 65.
- Collect information about lipid disorders in family members, including parents, siblings and offspring. There are three disorders to look for:

DYSLIPIDEMIA	Triglycerides greater than 150 mg/dL and HDL-C less than 40 mg/dL	
COMBINED HYPERLIPIDEMIA	Triglycerides greater than or equal to 150 mg/dL and LDL-C greater than 160 mg/dL, usually with HDL-C less than 40 mg/dL	
ISOLATED LOW HDL-C	Isolated HDL-C less than 40 mg/dL	

4

INTERPRET HDL MAP RESULTS

Initiate appropriate treatment after correcting contributing causes for low HDL-C and elevated triglycerides

PARAMETER	R/LAB VALUES	POTENTIAL DIAGNOSIS/CLINICAL SIGNIFICANCE	TREATMENT CONSIDERATIONS**	
Ť	ApoA-I in α -1 HDL>35 mg/dL	Optimal Indicates excellent reverse cholesterol transport	No treatment changes or additions for optimizing HDL particles	
Ŷ	ApoA-I in α -1 HDL>45 mg/dL			
	ApoA-I in α-1 HDL≥25 but≤35 mg/dL	Not optimal for CVD prevention	Lifestyle	
		· If CVD is not present, may consider treating only with	· Weight reduction if needed	
		lifestyle modification	· Exercise (>30 minutes/day)	
	ApoA-I in α-1 HDL≥35 but≤45 mg/dL	If CVD is present, consider treating to optimal levels	 Diet (low saturated fat (<7%), low cholesterol (<200 mg/day), low trans fat and low sugar) 	
		Absorbed by the second side of	· Smoking cessation	
M	ApoA-I in α -1 HDL<25 mg/dL	heart disease	Medications	
			 Use statin* to optimize LDL-C first; consider impact of statin on HDL particles^{4,7,8} 	
	ApoA-I in α -1 HDL<35 mg/dL		 Statin efficacy for beneficially modifying HDL in decreasing order: 	
ĥ	ApoA-I in α -1 HDL<25 mg/dL and	Abnormal levels associated with increased risk of beart disease	- Rosuvastatin · Atorvastatin · Simvastatin · Pravastatin · Lovastatin	
	prev- r ∕∠3 mg/u∟	Indicates a problem of HDL metabolism	- Niacin* (2 grams/day) is the best agent for optimizing HDL ^{1,9}	
Ĥ	ApoA-I in α-1 HDL<35 mg/dL and preß-1>25 mg/dL	Often associated with triglyceride levels >150 mg/dL	- Fenofibrate* is the best agent for lowering triglycerides, especially in those with TG >500 mg/dL; can also use fish oil (\geq 4 g/day).	

¹Nacin at a dose of 2 grams/day will increase apoA-I in large **a**-1 HDL by 100% or more. Potent statins such as rosuvastatin and atorvastatin will increase this parameter by 10-40%. Fibrates such as gemfibrozil or fenofibrate do not increase this parameter, but increase intermediate sized HDL particles such as **a** -3 and **a**-2 HDL. Statins are the best agents for lowering LDL and small dense LDL, naich is the best agent for traising large HDL, and fibrates are the best agents for triglyceride lowering. However all of the above agents will lower triglyceride levels. The statin/niacin combination has been shown to be the most effective tratammer regiment or promote regression of CDV or cardiovascut divasae. The space-1 level in very large **a**-1 HDL should be >20 mg/dL in men and >30 mg/dL in women.

IMPACT OF DRUG CLASSES ON HDL PARTICLES & METABOLISM

HDL PARAMETER	NIACIN ^{1,9}	STATINS 4, 7, 8	FIBRATES 10, 11
HDL-C	↑ 20–40%	↑ 2–10%	↑ 2–10%
ApoA-I concentration	1	_	-
α-1 particles	↑ up to 115%	↑ 12-36%	Slight ↓
α-3 particles	-	_	Slight ↑
Preß-1 particles	↓up to 30%	↓up to 40%	-
Metabolism	↑ apoA-I production ↑ ATP-binding cassette protein A1 (ABCA-1) expression in the liver	↓cholestryl ester transfer protein (CETP) activity No change in apoA-I kinetics	↑ gene expression of apoA-I, apoA-II & lipoprotein lipase (LPL) ↑ apoA-I fractional catabolic rate (FCR)

BOSTON HEART DIAGNOSTICS

Boston Heart Diagnostics is transforming the treatment of cardiovascular disease by providing healthcare providers and their patients with novel, personalized diagnostics and integrated customized lifestyle programs that have the power to change the way clinicians and patients communicate about disease and improve heart health. Boston Heart looks beyond the "good" and "bad" cholesterol assessment that conventional labs provide to give a more complete picture of heart health. Founded by renowned cardiovascular researchers and led by seasoned lab and diagnostic executives, Boston Heart is one of the fastest growing health companies in the country.

REFERENCES

- Asztalos BF, Batista M, Horvath KV, Cox CE, Dallal GE, Morse JS, Brown GB, Schaefer EJ. Change in α-1 HDL concentration predicts progression in coronary artery stenosis. Arterioscler Thromb Vasc Biol. 2003 May;23(5):847–852.
- Asztalos BF, Cupples LA, Demissie S, Horvath KV, Cox CE, Batista MC, Schaefer EJ. High-density lipoprotein subpopulation profile and coronary heart disease prevalence in male participants of the Framingham Offspring Study. Arterioscler Thromb Vasc Biol. 2004 Nov;24(11):2181–2187.
- Asztalos BF, Collins D, Cupples LA, Demissie S, Horvath KV, Bloomfield HE, Robins SJ, Schaefer EJ. Value of high-density lipoprotein (HDL) subpopulations in predicting recurrent cardiovascular events in the Veterans Affairs HDL Intervention Trial. Arterioscler Thromb Vasc Biol. 2005 Oct;25(10):2185–2191.
- Asztalos BF, LeMaulf F, Dallal GE, Stein E, Jones PH, Horvath KV, McTaggart F, Schaefer EJ. Comparison of the effects of high doses of rosuvastatin versus atorvastatin on the subpopulations of high-density lipoproteins. Am J Cardiol. 2007 Mar;99(5):681–685.
- Rubins HB, Robins SJ, Collins D, Nelson DB, Elam, MB, Schaefer EJ, Faas FH, Anderson JW; for the VA-HIT Study Group. Diabetes, plasma insulin, and cardiovascular disease: subgroup analysis from the Department of Veterans Affairs high-density lipoprotein intervention trial (VA-HIT). Arch Intern Med. 2002 Dec;162(22):2597–2604.
- Schaefer EJ, Lamon-Fava S, Ordovas JM, et al. Factors associated with low and elevated plasma high density lipoprotein cholesterol and apolipoprotein A-I levels in the Framingham Offspring Study. J. Lipid Res. 1994 May;35(5):871–882.

- Lamon-Fava S, Diffenderfer MR, Barrett PHR, Barrett HR, Buchsbaum A, Matthan NR, Lichtenstein AH, Dolnikowski GG, Horvath K, Asztalos BF, Zago V, Schaefer EJ. Effects of different doses of atorvastatin on human apolipoprotein B-100, B-48, and A-I metabolism. J Lipid Res. 2007 Aug;48(8):1746–1753.
- Asztalos BF, Horvath KV, McNamara JR, Roheim PS, Rubinstein JJ, Schaefer EJ. Comparing the effects of five different statins on the HDL subpopulation profiles of coronary heart disease patients. Atherosclerosis. 2002 Oct;164(2):361–369.
- Lamon-Fava S, Diffenderfer MR, Barrett PHR, Barrett HR, Buchsbaum A, Nyaku M, Horvath KV, Asztalos BF, Otokozawa S, Ai M, Matthan NR, Lichtenstein AH, Dolnikowski GG, Schaefer EJ. Extended-release niacin alters the metabolism of plasma apolipoprotein (Apo) A-I and ApoBcontaining lipoproteins. Arterioscler Thromb Vasc Biol. 2008 Sep;28(9):1672–1678.
- Watts GF, Barrett PH, Ji J, et al. Differential regulation of lipoprotein kinetics by atorvastatin and fenofibrate in subjects with the metabolic syndrome. Diabetes. 2003 Mar;52(3):803–811.
- Asztalos BF, Collins D, Horvath KV, Bloomfield HE, Robins SJ, Schaefer EJ. Relation of gemfibrozil treatment and high-density lipoprotein subpopulation profile with cardiovascular events in the Veterans Affairs High-Density Lipoprotein Intervention Trial. Metabolism. 2008 Jan;57(1):77–83.

877.425.1252 bostonheartdiagnostics.com

©2016 Boston Heart Diagnostics Corporation. All rights reserved. The Boston Heart logo and Boston Heart HDL Map are registered trademarks of Boston Heart Diagnostics Corporation. Boston Heart reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.

